SSPC1, SSPC2, SSPC3, SSPC4 Delivery Catheter

INSTRUCTIONS FOR USE - ENGLISH

Manufacturer
CenterPoint Systems LLC
3338 Parkway Blvd
West Valley City, UT 84119
USA
+1-877-848-0828
info@centerpoint-systems.com

Distributor
Boston Scientific Corp
4100 Hamline Ave N
St Paul, MN 55112
USA
+1-651-582-4000
boston.scientific.com

Device Description

The Delivery Catheter is a single-use percutaneous catheter intended to introduce various types of catheters and pacing or defibrillator leads. The Delivery Catheter is packaged with a dilator for introduction into the vasculature. Proximally, the Delivery Catheter is equipped with a low-friction hemostatic valve, and the distal tip, catheter shaft, and dilator are radiopaque to facilitate imaging under fluoroscopy. The Delivery Catheter has a braided shaft construction for increased torque response. The Delivery Catheter is designed to be slittable, thereby allowing its removal after device placement. A variety of curves are available to accommodate various anatomies and different lead locations. The Delivery Catheter has an inner diameter of 6.5F, an outer diameter of 8F, and the dilator is compatible with a 0.035" (0.89mm) guidewire. The catheter usable length is 40cm.

Indications for Use

The Delivery Catheter is indicated for the venous introduction of various types of catheters and pacing or defibrillator leads.

Contraindications

- Obstructed or inadequate vasculature for venous access

Warnings

- Do not advance Delivery Catheter against resistance without careful assessment of the cause of resistance under fluoroscopy.
- Do not advance a lead or catheter through the Delivery Catheter against resistance without careful assessment of the cause of resistance under fluoroscopy.
- Do not remove the Delivery Catheter or a lead or catheter placed through the Delivery Catheter against resistance without careful assessment of the cause of resistance under fluoroscopy.
- Do not resterilize or reuse. Structural integrity and/or function may be impaired by cleaning, resterilization, or reuse.
- Do not expose the device to the MR Environment

Precautions

- The Delivery Catheter should be used by interventionalists experienced in performing cardiovascular procedures.
- Care should be taken to consider vessel size and other relevant aspects of patient anatomy relative to the specified size and shape of the selected Delivery Catheter. Evaluations of this product were conducted using models representative of adult human anatomy.
- Do not use if package is open or damaged.
- Use by "Use by date."
- Exposure to temperatures above 54°C (130°F) may damage device and accessories.
- Upon removal from package, inspect device to ensure no damage has occurred during shipping.
- Do not expose to solvents.
- Use in conjunction with fluoroscopic guidance and proper anti-coagulation agents.
- Ensure the catheter is thoroughly flushed and free of air prior to use.
- Care should be taken during the procedure to ensure that the hemostatic valve is not damaged and that the side port remains closed to reduce potential for air ingress or blood loss.

Complications

Possible complications include, but are not limited to, the following: exposure to x-ray radiation, adverse or allergic reactions to contrast agents, infection, hematoma, pneumothorax, embolization, vessel thrombosis, dissection, acute occlusion, clot formation, hemorrhage, vessel rupture, arrhythmia or heart block, hemodynamic changes, myocardial infarction, perforation of the heart, cardiac tamponade, stroke, and death.

Device Preparation

- Flush the Delivery Catheter with saline.
- Place the dilator in the Delivery Catheter.
- Flush the dilator with saline.

Recommended Procedure

- Obtain and verify venous access.
- Advance the Delivery Catheter and dilator together over a guidewire and to the desired location using fluoroscopy.
- Remove the dilator and guidewire.

- Advance a lead or catheter through the Delivery Catheter to the desired location using fluoroscopy. Ensure that the outer diameter of the lead or catheter is less than the inner diameter of the Delivery Catheter before attempting to introduce the lead or catheter.
Note: For leads with a mannitol capsule, it may be necessary initiate the capsule dissolution process by placing the lead tip in the catheter hub and waiting for approximately 30 seconds until the lead is able to advance freely through the Delivery Catheter.
- When a lead or catheter has been positioned in the desired location, the Delivery Catheter may be removed by slitting if it is desired to leave the lead or catheter in place. Alternatively, the Delivery Catheter may be removed without slitting following standard vascular technique.
- To slit, stabilize the lead and slitter in one hand in a fixed position. With the other hand, retract the hub of the Delivery Catheter toward the slitter and cut through hub. Continue slowly retracting the Delivery Catheter from the vessel and in line with the slitter until the entire device has been slit and removed.

Device Disposal

- After use, dispose of the catheter and packaging in accordance with hospital, administrative and/or local government policy.

Symbols Glossary

	Medical device		Sterilized using radiation
	Contents		Do not store in direct sunlight
	Outer diameter		Keep dry
	Inner diameter		Do not use if package is damaged
	Maximum usable length		Do not reuse
	Maximum compatible guidewire outside diameter		Do not re-sterilize
	Sterile barrier system		Manufacturer
	Model number		Upper limit of temperature
	Lot number		Non-pyrogenic
	Date of manufacture		MR Unsafe
	Use by date		Distributor
	Consult Instructions for Use		EC authorized representative
	CAUTION: Federal (USA) law restricts this device for sale by or on the order of a physician.		

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY
THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE CENTERPOINT SYSTEMS PRODUCT(S) DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL CENTERPOINT SYSTEMS BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND CENTERPOINT SYSTEMS TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

Descriptions or specifications in CenterPoint Systems printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. CenterPoint Systems will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

EU Authorized Representative
ICON (LR) Limited
South County Business Park
Leopardstown Dublin 18 D18
X5R3 Ireland

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