CPS Locator™ 3D Delivery Catheter
INSTRUCTIONS FOR USE

Manufacturer
CenterPoint Systems LLC
3338 Parkway Blvd
West Valley City, UT 84119
USA
+1-877-848-0828
info@centerpoint-systems.com

Distributor
Abbott Medical
Cardiac Rhythm Management
Division
15900 Valley View Court
Sylmar, CA 91342 USA
+1 818 362 6822
800 722 3774 (US Only)

Device Description
The Delivery Catheter is a single-use percutaneous catheter intended to introduce pacing or defibrillator leads.

The Delivery Catheter is packaged with a dilator for introduction into the vasculature. Proximally, the Delivery Catheter is equipped with a hemostatic valve, and the distal tip is radiopaque to facilitate imaging under fluoroscopy. The Delivery Catheter is designed to be slittable (with slitters such as Agilis HisPro Slitter Model DS3A001, etc), thereby allowing its removal after device placement. A variety of curves (Small, Medium, Large) are available to accommodate various anatomies and different lead locations. The Delivery Catheter has an inner diameter of 7F (2.44mm), an outer diameter of 9F(3mm), and the dilator is compatible with a 0.035" (0.89mm) guidewire. The catheter comes in usable lengths of 42cm and 45cm.

Intended Purpose
The Delivery Catheter is intended for the venous introduction of pacing or defibrillator leads

Indications for Use
The Delivery Catheter is indicated for adult populations requiring pacemaker or defibrillator lead implantation for the treatment of heart arrhythmia. The Delivery Catheter is not indicated to treat or diagnose any disease or clinical condition, but is used to assist physicians achieve lead delivery via the right heart during lead implantation procedures.

Contraindications
• Obstructed or inadequate vasculature for venous access
• Review contraindications associated with ancillary devices used with the Delivery Catheter

Warnings
• Do not advance Delivery Catheter against resistance without careful assessment of the cause of resistance under fluoroscopy.
• Do not advance a lead through the Delivery Catheter against resistance without careful assessment of the cause of resistance under fluoroscopy.
• Do not remove the Delivery Catheter or a lead placed through the Delivery Catheter against resistance without careful assessment of the cause of resistance under fluoroscopy.
• Do not resterilize or reuse. Structural integrity and/or function may be impaired by cleaning, resterilization, or reuse.
• Do not expose the device to the MR Environment

Precautions
• The Delivery Catheter should be used by interventionalists experienced in performing cardiovascular procedures.
• Care should be taken to consider vessel size and other relevant aspects of patient anatomy relative to the specified size and shape of the selected Delivery Catheter. Evaluations of this product were conducted using models representative of adult human anatomy.
• Do not use if package is open or damaged.
• Keep stored in shelf carton until use.
• Use by "Use by date."
• Store at room temperature. Data from the accelerated aging support a recommended continuous storage temperature of 23°C. Based on the same data, the product may withstand transient excursions of up to 54°C.
• Upon removal from package, inspect device to ensure no damage has occurred during shipping.
• Do not expose to solvents.
• Use in conjunction with fluoroscopic guidance and proper anti-coagulation agents.
• Ensure the catheter is thoroughly flushed and free of air prior to use.
• Care should be taken during the procedure to ensure that the hemostatic valve is not damaged and that the side port remains closed to reduce potential for air ingress or blood loss.

Complications
Possible complications associated with the Delivery Catheter include: delay in procedure, vessel trauma, endocardial trauma, air ingress, blood loss, lead dislodgement, embolization (material, air), allergic reaction, lead dislodgement, infection, device incompatibility, incomplete therapy.
Possible complications associated with the use of ancillary devices with the Delivery Catheter include: Air embolism, Bleeding, Bradycardia, Breakage/failure of the implant instruments, Cardiac perforation, Cardiac tamponade, Chronic nerve damage, Death, Extracardiac stimulation (muscle/nerve stimulation), Hematoma, Heart block, Hemothorax, Incisional pain, Infection including endocarditis, Lead dislodgement, Myocardial infarction, Myocardial trauma, Pericardial effusion, Pericarditis, Pleural Effusion, Pneumothorax, Stroke, Tachyarrhythmias, Thrombosis/thromboemboli, Valve damage, Venous occlusion, Venous trauma (e.g., perforation, dissection, erosion)

Patient Target Group
Adult populations requiring pacemaker/defibrillator implantation
Intended users
Interventionalists experienced in performing cardiovascular procedures and who have access to the appropriate facilities (fluoroscopy) to perform the procedures

Clinical Benefits, Performance Characteristics
The Delivery Catheters are intended to be used during pacemaker/defibrillator lead implantation procedures to assist with lead implant. By helping to achieve optimal lead placement during implantation the intended physiologic effect of pacing is more likely to be achieved. The overall clinical benefits of pacing therapy include management of cardiac arrhythmias, a decrease in symptoms (e.g. syncope, dizziness, fatigue, shortness of breath, chest pains), decreased dependence on medications, reduced cost of care, increased exercise capacity, and an overall increase in quality of life.

Summary of Safety and Clinical Performance
For customers in the European Union, use the device name found in the labeling to search for the device's Summary of Safety and Clinical Performance, which is available on the European database on medical devices (Eudamed) website:
https://ec.europa.eu/tools/eudamed

Until Eudamed becomes operational, the Summary of Safety and Clinical Performance (SSCP) can be obtained upon request to the manufacturer. To request a copy of the SSCP, contact CenterPoint Systems at info@centerpoint-systems.com.

Device Preparation
1. Flush the Delivery Catheter with saline.
2. Place the dilator in the Delivery Catheter.
3. Flush the dilator with saline.

Recommended Procedure
1. Obtain and verify venous access.
2. Advance the Delivery Catheter and dilator together over a guidewire and to the desired location using fluoroscopy.
3. Remove the dilator and guidewire.
4. Advance a lead through the Delivery Catheter to the desired location using fluoroscopy. Ensure that the outer diameter of the lead is less than the inner diameter of the Delivery Catheter before attempting to introduce the lead.
5. When a lead has been positioned in the desired location, the Delivery Catheter may be removed by slitting if it is desired to leave the lead in place. Alternatively, the Delivery Catheter may be removed without slitting following standard vascular technique.
6. To slit, stabilize the lead and slitter in one hand in a fixed position. With the other hand, retract the hub of the Delivery Catheter toward the slitter and cut through hub. Continue slowly retracting the Delivery Catheter from the vessel and in line with the slitter until the entire device has been slit and removed.

Device Disposal
It is important to understand and follow all laws regarding the safe and proper disposal of medical devices as infection or other biological harm could occur. Dispose according to standard operational procedure for biohazardous/potential sharps-containing products.

Complaint Reporting
Any complaints and/or serious incidents involving CenterPoint Systems products should be reported to CenterPoint Systems LLC and the competent authority of the Member State in which the user and/or patient is established immediately. To report a complaint, contact CenterPoint Systems at complaints@centerpoint-systems.com.

Table with 2 columns: Symbol/Icon and Description. Includes symbols for Medical device, Contents, Inner Diameter/Outer diameter, Unique Device Identifier, Sterilized using radiation, Do not store in direct sunlight, Keep dry, Do not use if package is damaged, Do not reuse, Do not re-sterilize, Manufacturer, Upper limit of temperature, Non-pyrogenic, MR Unsafe, EC authorized representative, Lot number, Date of manufacture, Use by date, Consult Instructions for Use, One device included in the package, Translation performed by Lionbridge Technologies, LLC.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY
THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE CENTERPOINT SYSTEMS PRODUCT(S) DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL CENTERPOINT SYSTEMS BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND CENTERPOINT SYSTEMS TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN. Descriptions or specifications in CenterPoint Systems printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. CenterPoint Systems will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

EU Authorized Representative
MedEnvoy Global BV
Prinses Margrietplantsoen 33
Suite 123
2595 AM Den Haag
Netherlands
SRN: NL-AR-000024028

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Çeviriyi gerçekleştiren taraf: Lionbridge Technologies, LLC 890 Winter Street, Waltham, MA 02451
AB Yetkili Temsilcisi