INSTRUCTIONS FOR USE – ENGLISH

Helix Locking Tool (5616)

Description

Package contains two (2) Helix Locking Tools.

Tool is compatible with IS-1¹ lead connectors and lead stylets up to 0.016" (0.4mm) diameter.

The Helix Locking Tool is a sterile, single-use device compatible with pacing and/or defibrillator leads with an extendable-retractable fixation helix and an IS-1 connector with terminal pin. The Helix Locking Tool is press-fit to the lead terminal pin to fix in place or, upon rotation, transfer torque to the terminal pin. The lead helix may be extended or retracted based on the direction of torque applied. The tool may also guide stylet insertion into the lead.

For additional information regarding appropriate use of accessories with a particular lead (including any applicable warnings, precautions, and other safety information), consult the literature provided by the lead manufacturer.

Intended Use/Indication for Use

The Helix Locking Tool is intended to facilitate lead fixation by providing control over the extension and retraction of the lead helix of an IS-1 pacing or defibrillator lead. It is also intended to facilitate insertion of a stylet into the lead through the terminal pin.

Warnings

- When using the tool, do not rotate the lead terminal pin more than the maximum number of turns recommended by the lead manufacturer. Continuing to rotate the terminal pin can cause lead damage, lead dislodgement, tissue trauma, and/or voltage thresholds to rise.
- Remove the tool from the lead terminal pin before placing the lead terminal pin in the patient and securing the lead and/or connecting the lead to a pacemaker/defibrillator.
- For single use/single procedure only. Do not reuse, reprocess, or. Such may cause failure of the tool and/or infection or cross-contamination.
- Do not expose the tool to the MR environment. Testing has not been performed to establish safe use in this environment.

Precautions

- Read these instructions and consult the literature provided by the lead manufacturer prior to use.
- Inspect the tool and its packaging prior to use. Do not use if the package is open or there is any evidence of package or tool damage.
- Store in a cool, dark, dry location.
- Store at room temperature. Data from accelerated aging supports a recommended continuous storage temperature of 23°C. Based on the same data, the product may withstand transient temperature excursions of up to 54°C.

Complications

Possible complications associated with use of the Helix Locking Tool include, but are not limited to, delay in procedure, lead dislodgement, device incompatibility and/or incomplete therapy.

Patient Target Group


Adult populations requiring pacing or defibrillator implantation.

Intended Users

Medical professionals experienced in performing pacing and/or defibrillator implant procedures who have access to the appropriate facilities to perform such procedures.

Directions

- Press-fit the Helix Locking Tool to the terminal pin of a lead with an IS-1 connector. The arrow on the tool body should point toward the lead.

- Rotate the tool (and lead terminal pin) clockwise to extend or counterclockwise to retract the lead helix to facilitate lead fixation. Do not exceed the maximum number of turns recommended by the lead manufacturer.
- A stylet having diameter 0.016" (0.4mm) or less may pass through the tool into the lead terminal pin.
- To remove the tool, pull it gently away from the lead terminal pin.

Disposal

All items used during invasive procedures, including the tool and its packaging, may be contaminated with potentially infectious material. **CAUTION:** Dispose using proper biohazard handling technique in accordance with institutional policy and/or applicable law.

Complaint Reporting

If you have reason to believe that a serious incident occurred while using this device, please report it to the manufacturer. In the European Union, report this also to your national authority.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THE PRODUCTS DESCRIBED HEREIN ARE PROVIDED "AS IS" WITH ALL FAULTS. THERE IS NO EXPRESS OR IMPLIED WARRANTY FOR SUCH PRODUCTS INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. UNDER NO CIRCUMSTANCES SHALL CENTERPOINT SYSTEMS BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND CENTERPOINT SYSTEMS TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

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